Collagen-Based Bone Filling Augmentation Material

Instructions for Use

Please read these instructions carefully before using FOUNDATION.

Product Overview

FOUNDATION is made of absorbable atelo-collagen. The atelo-collagen has been cross-linked by heat treatment in order to achieve biocompatibility. FOUNDATION is provided in bullet shapes. The collagen (liquid form) suspension is pored into a bullet shape mold. The design of the bullet shape is to permit easy placement into the extraction socket. It consists of fibrillar and heat-denatured atelo-collagen. Fibrillar atelo-collagen provides scaffolding for surrounding cells and the heat-denatured atelo-collagen stimulates infiltration of the cells in the product. Contents are sterile and non-pyrogenic.

Indication for Use

FOUNDATION is a collagen-based bone filling augmentation material for use in the filling of extraction sockets.

Contraindication

This product uses proteins of bovine origin as raw material and should not be used in those patients who have a known history of hypersensitivity.

Warnings

1. Use of this product should be discontinued immediately if the patient displays any symptoms such as rash, hives, or other allergic reactions.
2. Do not touch the material with your bare hands.

Precautions

1. Do not resterilize or reuse. This product is sterile and intended for single use only.
2. Do not use if the package and/or product is damaged.
3. Once the unit package has been opened, discard if not used immediately.
4. Use of this product in patients predisposed to asthma, hives, or other allergic reactions should be evaluated prior to use.
5. Use of this product in patients who have just given birth or who are nursing a baby should be evaluated prior to use.
6. The safety of this product in pregnant patients has not been studied or established, therefore, use in pregnant patients should be evaluated for associated risks.
7. This product has no antimicrobial properties. Any decision regarding the use of an antibiotic or other therapies when using this product should be made based upon the patient's risk of infection.
8. Applying excessive amounts of the material may cause tenderness and irritation of the socket.
9. Secure FOUNDATION by suturing, or by some other means. Failure to do so can cause dislodgement of the material which can present a choking hazard.

Caution: US Federal Law restricts this device to sale by or on the order of a dentist.

How to Open the Container

1. Open aluminum package at TEAR HERE and remove plastic container.
2. While holding the edge of the container, remove FOUNDATION using a sterile instrument. Do not touch FOUNDATION with your bare hands.

Instructions for Use

1. If the extraction socket is bleeding excessively, achieve hemostasis with gauze or other appropriate material.
2. Curette the socket to expose bone, and if needed, curette or abrade the surrounding gingiva to stimulate bleeding to ensure there will be infiltration of blood into the product.
 - a) If the product is to be used in an extraction socket with an oroantral fistula that has been left untreated for an extended period of time resulting in the oral mucous epithelium covering the bone surface of the extraction socket with the fistula, strip the oral mucous epithelium from the bone surface of the extraction socket with the fistula before insertion of FOUNDATION.
 - b) If the product is to be used in a socket with chronic inflammation or other disease, curette the socket to remove any granulation tissue before application of FOUNDATION. Failure to do so may dislodge FOUNDATION from the socket.
3. Select a slightly larger size of FOUNDATION than the socket to be filled to ensure direct contact with bone. Remove FOUNDATION from the container and place it in the socket with a sterile instrument.
4. Secure FOUNDATION by suturing, or some other means, to prevent its dislodgement from the socket.
5. As necessary, immediate rinsing, gargling, or other appropriate treatment should be considered to minimize the risk of infection.

Storage

Keep dry. Keep away from heat. Avoid direct sunlight. Do not store under conditions of extreme humidity or room temperatures exceeding 86° F (30° C).

Packaging

Item/Product Number	Dimensions (dia. x length)	No. of packages (per box)
S (Small) 27-500-100	Ø 8 x 25 mm	10
M (Medium) 27-500-200	Ø 15 x 25 mm	5
Assortment Pack 27-500-150	(includes both)	3 Small, 3 Medium

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